



Podology abrasive caps and abrasive cap carriers




New shapes for even higher precision!

We have developed new shapes for our innovative grinding caps and grinding cap carriers, to enable us to provide even more suitable tools for podiatric applications. As well our patented SK THERMO as our SK PODO are now available in the „conical“ and „pointed“ versions in the qualities coarse, medium and fine. The matching grinding cap carriers are available in different sizes and also in „conical“ and „pointed“ finish. Our new shapes are CE marked and meet the legal requirements for medical products.








SK THERMO – Extra gentle protection with color change technology

Our state-of-the-art THERMO abrasive caps warn against overheating and damage to the skin through a 2-step color change. So they provide more safety in medical foot care and a pleasant treatment without tingling and pain.

SK THERMO, conical

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A2108112080	SK THERMO 11 mm conical	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A2108112150	SK THERMO 11 mm conical	medium	7.000 – 10.000 (max. 11.000)	5 x 10







SK THERMO, pointed

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21080073080	SK THERMO 7 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080073150	SK THERMO 7 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080103080	SK THERMO 10 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080103150	SK THERMO 10 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080133080	SK THERMO 13 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21080133150	SK THERMO 13 mm spitz pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10












SK PODO – Hygienic protection with non-stick effect

The SK PODO can be used to remove large quantities of callused skin even with low contact pressure. Thanks to its innovative design and non-stick effect, the sanding cap is significantly gentler and more efficient to use, since practically no skin residues are deposited on the cap surface.

SK PODO, conical

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A2110051080	SK PODO 5 mm conical	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A2110051150	SK PODO 5 mm conical	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A2110051320	SK PODO 5 mm conical	fine	7.000 – 10.000 (max. 11.000)	5 x 10
		A2110112080	SK PODO 11 mm conical	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A2110112150	SK PODO 11 mm conical	medium	7.000 – 10.000 (max. 11.000)	5 x 10

SK PODO, pointed

Product picture	1:1	Product number	Description	Quality	Recommended speed/rpm	
		A21100053080	SK PODO 5 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100053150	SK PODO 5 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100073080	SK PODO 7 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100073150	SK PODO 7 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100103080	SK PODO 10 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100103150	SK PODO 10 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100133080	SK PODO 13 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100133150	SK PODO 13 mm pointed	medium	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100163080	SK PODO 16 mm pointed	coarse	7.000 – 10.000 (max. 11.000)	5 x 10
		A21100163320	SK PODO 16 mm pointed	fine	7.000 – 10.000 (max. 11.000)	5 x 10









* Every box contains 5 PE bags with 10 abrasive caps each.



Abrasive cap carriers GT PODO – precise and custom-fit

The GT PODO grinding cap carrier guarantees high accuracy of fit and safe handling even in the pointed and conical shape.

GT PODO, Shank 2,35 x 40 mm

Product picture	1:1	Product number	Description	Recommended speed/rpm	 *
		A25100512	GT PODO 5 mm conical	7.000–10.000 (max. 11.000)	5 x 1
		A25101122	GT PODO 11 mm conical	7.000–10.000 (max. 11.000)	5 x 1
		A251000532	GT PODO 5 mm pointed	7.000–10.000 (max. 11.000)	5 x 1
		A251000732	GT PODO 7 mm pointed	7.000–10.000 (max. 11.000)	5 x 1
		A251001032	GT PODO 10 mm pointed	7.000–10.000 (max. 11.000)	5 x 1
		A251001332	GT PODO 13 mm pointed	7.000–10.000 (max. 11.000)	5 x 1
		A251001632	GT PODO 16 mm pointed	7.000–10.000 (max. 11.000)	5 x 1

* Every box contains 5 PE bags with 1 abrasive cap carrier each.

User-friendly, customer-oriented, innovative and flexible:
LUKAS Podological – Specialist for abrasive caps and carriers

LUKAS-ERZETT Vereinigte Schleif- und
 Fräswerkzeugfabriken GmbH & Co. KG
 Gebrüder-Lukas-Straße 1
 51766 Engelskirchen (Germany)
 Fon +49 2263 84-0
 Fax +49 2263 84-327 (Domestic Sales)
 Fax +49 2263 84-300 (International Sales)
 le@lukas-erzett.de

Your local reseller: